

A Comprehensive Guide to Time Domain Feature Extraction and Time-Frequency Analysis

Time Domain Feature Extraction (Step 3)

This phase focuses on extracting relevant features directly from a signal's amplitude over time, requiring careful steps for accurate results.

Input Data

The data used for time domain feature extraction can be either **raw** or **filtered**. The choice between these depends on the earlier processing steps to enhance the signal for better analysis.

Windowing Techniques



Applying windowing techniques is a fundamental step in segmenting data effectively. The MLS offers two main types of windowing:


1. Non-Overlapping Windowing:

- Windows are applied separately without any overlap.
- Ideal for data that can be divided into independent segments or for signals with abrupt changes.
- Suitable when continuous analysis isn't required.

2. Overlapping Windowing:

- Windows overlap to provide a more detailed and continuous analysis.
- Particularly suited for continuous signals or when gradual changes are significant.
- Commonly used in bio-signals like EEG, ECG, EMG, and in speech analysis.





Best Practice

The choice between non-overlapping and overlapping windowing hinges on the data type and the analysis objective. For **complex signals**, **overlapping windowing** with an appropriate window size is recommended. This should be adjusted based on the *sampling rate* and signal characteristics for higher accuracy.

Common Feature Methods

The software includes **22 different methods** for feature extraction. Some key methods are:

- **Mean:** Provides the average value of the signal within a window.
- **Absolute Standard Deviation:** Measures the variability of the signal.
- **Root Mean Square (RMS):** Calculates the square root of the mean of the squares of the signal values, useful for assessing the power of the signal.

Management & Visualization

Once a feature is selected, it can be managed by saving or deleting. The "Display Plot" option allows for visual analysis of the extracted features, which aids in interpretation and understanding of the data.

Saving Features

To facilitate direct use in classification tasks, extracted features should be saved using the "Save Feature" option.

Time-Frequency Domain Analysis: Wavelet Transformation (Step 4)

Wavelet Transformation is essential for analyzing complex, especially non-stationary signals across both time and frequency domains.

Suitability

This method is particularly powerful for **non-stationary signals** such as bio-signals (ECG, EEG) and audio signals.

Input Data

For wavelet transformation, you can upload raw data, filtered data, or even feature signals.

Core Purposes

Wavelet transformation serves multiple purposes:


- **Data Compression:** Reduces noise while retaining key information.
- **Multi-Resolution Analysis:** Examines the signal at different levels of detail, separating it into Approximation and Detail components.
- **Noise Removal (De-noising):** Achieved using specific threshold rules (e.g., **rigsure**) and threshold types (e.g., **s**).

Key Parameters for Quality

- **Wavelet Type:** The choice of wavelet impacts analysis quality. The "haar" wavelet is simple and fast, suitable for discrete signals with low detail. For more complex signals, "db4" (Daubechies 4) or "sym4" (Symlet 4) wavelets are recommended as they preserve more details.
- **Decomposition Level (Level):** Indicates the depth of signal breakdown. Level 1 is sufficient for simple signals, while levels 3 to 5 can provide more accurate information for complex signals with high noise. However, higher levels increase computation time.

Visualization & Saving

The "Display Plot" option allows visualization of the decomposition results and interpretation of wavelet coefficients. After processing, wavelet coefficients can be saved using the "Save Coefficients" option for further use. These coefficients provide powerful time-frequency features for subsequent analysis.



🔗 Explore the repository here: <https://lnkd.in/dDJUyNtc>

Stay updated by following @rezasaadatyar